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REFERENCIAS PARTE 1

- Brennan N. Corneal oxygenation during contact lens wear: comparison of diffusion and EDP-based flux models. *Clin Exp Optom* 2005; 88:103-108.
- Brennan N. Beyond flux: total corneal oxygen consumption as an index of corneal oxygenation during contact lens wear. *Optom Vis Sci* 2005; 82:467-472.
- Brennan NA. A model of oxygen flux through contact lenses. *Cornea* 2001; 20:104-108.
- Compan V, Lopez-Alemay A, Riande E, et al. Biological oxygen apparent transmissibility of hydrogel contact lenses with and without organosilicon moieties. *Biomater* 2004; 25:359-365.
- Huang P, Zwang-Weissman J, Weissman BA. Is contact lens "T" still important? *Contact Lens Ant Eye* 2004; 27:9-14.
- Radke CJ, Chhabra M. Minimum Contact Lens Oxygen Transmissibility (Dk/L) With Monod Kinetics for the Corneal Oxygen Consumption Rate. *Invest. Ophthalmol. Vis. Sci.* 2005; 46:904-.
- Fonn D, Bruce A. A review of the Holden-Mertz criteria for critical oxygen transmission. *Eye Contact Lens* 2005; 31:247-251.
- Fonn D, Sweeney D, Holden B, et al. Corneal oxygen deficiency. *Eye Contact Lens* 2005; 31:23-27.
- Holden BA, Stretton S, Lazon de la Jara P, et al. The future of contact lenses: Dk really matters. *Contact Lens Spectrum* 2006; 21:20-28.
- Bonanno JA, Clark C, Pruitt J, et al. Tear oxygen tension under silicone-hydrogel lenses in human subjects. *Invest Ophthalmol Vis Sci* 2008; 49:E-abstract 4846.
- Fatt I. New physiological paradigms to assess the effect of lens oxygen transmissibility on corneal health. *CLAOJ* 1996; 22:25-29.
- Fatt I, Freeman RD, Lin D. Oxygen tension distributions in the cornea: a re-examination. *Exp Eye Res* 1974; 18:357-365.
- Bruce AS. Local oxygen transmissibility of disposable contact lenses. *Contact Lens Anterior Eye* 2003; 26:189-196.
- Efron N, Morgan PB, Cameron ID, et al. Oxygen permeability and water content of silicone hydrogel contact lens materials. *Optom Vis Sci* 2007; 84:328-337.
- Brennan N. Corneal oxygen consumption beneath toric soft lenses. 2008: In preparation.
- Leibowitz H, Laing R, Sandstrom M. Continuous wear of hydrophilic contact lenses. *Arch Ophthalmol* 1973; 89:306-310.
- Sarver MD, Staroba JE. Corneal edema with contact lenses under closed-eye conditions. *Am J Optom Physiol Opt* 1978; 55:739-743.
- Binder P. Complications associated with extended wear of soft contact lenses. *Ophthalmol* 1979; 86:1093-1101.
- Holden B, Sweeney D, Swarbrick H, et al. The vascular response to long-term extended contact lens wear. *Clin Exp Optom* 1986; 69:112-119.
- Zantos S, Holden B. Transient endothelial changes soon after wearing soft contact lenses. *Am J Optom Physiol Opt* 1977; 54:856-858.
- Zantos S, Holden B. Ocular changes associated with continuous wear of contact lenses. *Aust J Optom* 1978; 61:418-426.
- Zantos S. Cystic formations in the corneal epithelium during extended wear of contact lenses. *Int Contact Lens Clin* 1983; 10:128-146.
- Holden B, Sweeney D, Vannas A, et al. Effects of long-term extended contact lens wear on the human cornea. *Invest Ophthalmol Vis Sci* 1985; 26:1489-1501.
- Schoessler J. Corneal endothelial polymegathism associated with extended wear. *Int Contact Lens Clin* 1983; 10:148-156.
- Dumbleton KA, Chalmers RL, Richter DB, et al. Changes in myopic refractive error with nine months' extended wear of hydrogel lenses with high and low oxygen permeability. *Optom Vis Sci* 1999; 76:845-849.
- Rubén M, Brown N, Lobascher D, et al. Clinical manifestations secondary to soft contact lens wear. *Br J Ophthalmol* 1976; 60:529-531.
- McMonnies CW, Chapman-Davies A, Holden BA. The vascular response to contact lens wear. *Am J Optom Physiol Opt* 1982; 59:795-799.
- MacRae S, Matsuda M, Shelland S, et al. The effects of hard and soft contact lenses on the endothelium. *Am J Ophthalmol* 1986; 102:50-57.
- Farris R, Donn A. Corneal respiration with soft contact lenses. *J Am Optom Assoc* 1972; 43:292-294.
- Roth HW. The etiology of ocular irritation in soft lens wearers: distribution in a large clinical sample. *Contact Lens Intraoc Lens Med J* 1978; 4:38.
- McMonnies C. Contact lens-induced corneal vascularization. *Int Contact Lens Clin* 1983; 10:12-21.
- Efron N. Vascular response of the cornea to contact lens wear. *J Am Optom Assoc* 1987; 58:836-846.
- Holden B, Sweeney D, Sanderson G. The minimum precorneal oxygen tension to avoid corneal oedema. *Invest Ophthalmol Vis Sci* 1984; 25:476-480.
- Polse K, Mandell R. Critical oxygen tension at the corneal surface. *Arch Ophthalmol* 1971; 84:505-508.
- Holden B, Mertz G. Critical oxygen levels to avoid corneal edema for daily and extended wear contact lenses. *Invest Ophthalmol Vis Sci* 1984; 25:1161-1167.
- La Hood D, Grant T. Striae and folds as indicators of corneal oedema. *Optom Vis Sci* 1990; 67 (suppl.):196.
- Morgan PB, Maldonado-Codina C, Quhill W, et al. Central and peripheral oxygen transmissibility thresholds to avoid corneal swelling during open eye soft contact lens wear. *J Biomed Mat Res* 2008: Submitted for publication.
- Brennan NA, Fonn D. Corneal hypoxia, in *Siliconehydrogels: continuous wear contact lenses*, D.F. Sweeney, Editor. 2004, Butterworth-Heinemann: London.
- Cox I, Zantos S, Orsborn G. The overnight corneal swelling response of non-wear, daily wear, and extended wear soft lens patients. *Int Contact Lens Clin* 1990; 17:134-137.
- Brennan NA, Coles M-LC. Continuous wear, in *Contact Lens Practice*, N. Efron, Editor. 2002, Butterworth-Heinemann: Oxford. p. 275-294.
- Fonn D, du Toit R, Simpson TL, et al. Sympathetic swelling response of the control eye to soft lenses in the other eye. *Invest Ophthalmol Vis Sci* 1999; 40:3116-3121.
- Holden BA, Grant T, Kotow M, et al. Epithelial microcysts with daily and extended wear of hydrogel and rigid gas permeable contact lenses. *Invest Ophthalmol Vis Sci* 1987; 28(Suppl):372.
- Hickson S, Papas E. Prevalence of idiopathic corneal anomalies in a non contact lens-wearing population. *Optom Vis Sci* 1997; 74:293-297.
- Holden BA, Sweeney DF. The significance of the microcyst response: a review. *Optom Vis Sci* 1991; 68:703-707.
- Sweeney DF, Keay L, Jalbert I, et al. Clinical performance of silicone hydrogel lenses, in *Silicone hydrogels: the rebirth of continuous wear contact lenses*, D.F. Sweeney, Editor. 2000, Butterworth-Heinemann: Oxford. p. 90-149.
- Brennan NA, Coles ML, Connor HR, et al. A 12-month prospective clinical trial of comfilcon A silicone-hydrogel contact lenses worn on a 30-day continuous wear basis. *Contact Lens Anterior Eye* 2007; 30:108-118.
- Dumbleton KA, Chalmers RL, Richter DB, et al. Vascular response to extended wear of hydrogel lenses with high and low oxygen permeability. *Optom Vis Sci* 2001; 78:147-151.
- Brennan NA, Coles M-LC, Levy B, et al. One-Year Prospective Clinical Trial of balafilcon A (PureVision) Silicone-hydrogel Contact Lenses Used on a 30-Day Continuous Wear Schedule. *Ophthalmol* 2002; 109:1172-1177.

REFERENCIAS PARTE 1

- Santodomingo-Rubido J, Wolffsohn JS, Gilmartin B. Changes in ocular physiology, tear film characteristics, and symptomatology with 18 months silicone hydrogel contact lens wear. *Optom Vis Sci* 2006; 83:73-81.
- Santodomingo-Rubido J, Wolffsohn JS, Gilmartin B. Adverse events and discontinuations during 18 months of silicone hydrogel contact lens wear. *Eye Contact Lens* 2007; 33:288-292.
- Papas E. On the relationship between soft contact lens oxygen transmissibility and induced limbal hyperaemia. *Exp Eye Res* 1998; 67:125-131.
- Brennan NA, Coles ML, Ang JH. An evaluation of silicone-hydrogel lenses worn on a daily wear basis. *Clin Exp Optom* 2006; 89:18-25.
- Guillon M, Maissa C. The effect of optimising silicone hydrogels key physical properties on extended wear clinical performance. *Invest Ophthalmol Vis Sci* 2006; 49:ARVO E-abstract 2382.
- Maldonado-Codina C, Morgan PB, Schnider CM, et al. Short-term physiological response in neophyte subjects fitted with hydrogel and silicone-hydrogel contact lenses. *Optom Vis Sci* 2004; 81:911-921.
- Morgan PB, Efron N. Comparative clinical performance of two silicone hydrogel contact lenses for continuous wear. *Clin Exp Optom* 2002; 85:183-192.
- Papas E, Willcox MDP. Reducing the consequences of hypoxia: the ocular redness response. *Contact Lens Spectrum* 2006; 21 (suppl.):32-37.
- Inagaki Y, Akahori A, Sugimoto K, et al. Comparison of corneal endothelial bleb formation and disappearance processes between rigid gas-permeable and soft contact lenses in three classes of dk/l. *Eye Contact Lens* 2003; 29:234-237.
- Brennan NA, Coles M-LC, Connor RM, et al. Short-term corneal endothelial response to wear of silicone hydrogel contact lenses in east asian eyes. *Eye Contact Lens* 2008; 34:317-321.
- Chang SW, Hu FR, Lin LL. Effects of contact lenses on corneal endothelium -a morphological and functional study. *Ophthalmologica* 2001; 215:197-203.
- Yamauchi K, Hirst LW, Enger C, et al. Specular microscopy of hard contact lens wearers II. *Ophthalmology* 1989; 96:1176-1179.
- Carlson KH, Bourne WM. Endothelial morphologic features and function after long-term extended wear of contact lenses. *Arch Ophthalmol* 1988; 106:1677-1679.
- Stocker EG, Schoessler JP. Corneal endothelial polymegathism induced by PMMA contact lens wear. *Invest Ophthalmol Vis Sci* 1985; 26:857-863.
- Schoessler JP, Barr JT, Freson DR. Corneal endothelial observations of silicone elastomer contact lens wearers. *Int Contact Lens Clin* 1984; 11:337.
- Liu Z, Pflugfelder SC. The effects of long-term contact lens wear on corneal thickness, curvature, and surface regularity. *Ophthalmology* 2000; 107:105-111.
- Myrowitz EH, Melia M, O'Brien TP. The relationship between long-term contact lens wear and corneal thickness. *Clao J* 2002; 28:217-220.
- Alharbi A, Swarbrick HA. The effects of overnight orthokeratology lens wear on corneal thickness. *Invest Ophthalmol Vis Sci* 2003; 44:2518-2523.
- Binder P. Myopic extended wear with the Hydrocurve II soft contact lens. *Ophthalmol* 1983; 90:623-626.
- Harris MG, Sarver MD, Polse KA. Corneal curvature and refractive error changes associated with wearing hydrogel contact lenses. *Am J Optom Physiol Opt* 1975; 52:313-319.
- Fonn D, MacDonald KE, Richter D, et al. The ocular response to extended wear of a high Dk silicone hydrogel contact lens. *Clin Exp Optom* 2002; 85:176-182.
- Patel S. Changes in myopic refractive error with nine months' extended wear of hydrogel lenses with high and low oxygen permeability. *Optom Vis Sci* 2000; 77:285.
- Smith EL, Greeman P, Ho A, et al. Methods and apparatuses for altering relative curvature of field and positions of peripheral off-axis focal positions. US patent 7025460 2006.
- Szczotka-Flynn LB. Unintended Ortho-k Effects From Silicone Hydrogel Lenses. *Contact Lens Spectrum* 2004; 19(8):19-20.

REFERENCIAS PARTE 2

- Morgan PB, Efron N, Hill EA, et al. Incidence of keratitis of varying severity among contact lens wearers. *Br J Ophthalmol* 2005; 89:430-436.
- Schein OD, McNally JJ, Katz J, et al. The incidence of microbial keratitis among wearers of a 30-day silicone hydrogel extended-wear contact lens. *Ophthalmology* 2005; 112:2172-2177.
- Stapleton F, Keay L, Edwards K, et al. The incidence of contact lens-related microbial keratitis in Australia. *Ophthalmology* 2008; 115:1655-1662.
- Dart JK, Radford CF, Minassian D, et al. Risk factors for microbial keratitis with contemporary contact lenses: a case-control study. *Ophthalmology* 2008; 115:1647-1654.
- Morgan PB, Efron N, Brennan NA, et al. Risk factors for the development of corneal infiltrative events associated with contact lens wear. *Invest Ophthalmol Vis Sci* 2005; 46:31363143.
- Stapleton F, Keay L, Sanfilippo PG, et al. Relationship between climate, disease severity, and causative organism for contact lens-associated microbial keratitis in Australia. *Am J Ophthalmol* 2007; 144:690-698.
- Keay L, Edwards K, Naduvilath T, et al. Factors affecting the morbidity of contact lens-related microbial keratitis: a population study. *Invest Ophthalmol Vis Sci* 2006; 47:4302-4308.
- Szczotka-Flynn L, Diaz M. Risk of corneal inflammatory events with silicone hydrogel and low dk hydrogel extended contact lens wear: a meta-analysis. *Optom Vis Sci* 2007; 84:247-256.
- Andrasko G, Ryan K. Corneal staining and comfort observed with traditional and silicone hydrogel lenses and multipurpose solution combinations. *Optometry* 2008; 79:444-454.
- O'Leary DJ, Wilson G, Henson DB. The effect of anoxia on the human corneal epithelium. *Am J Optom Physiol Opt* 1981; 58:472-476.
- Skotnitsky CC, Naduvilath TJ, Sweeney DF, et al. Two presentations of contact lens-induced papillary conjunctivitis (CLPC) in hydrogel lens wear: local and general. *Optom Vis Sci* 2006; 83:27-36.
- Spring TF. Reaction to hydrophilic lenses. *Med J Aust* 1974; 1:449-450.
- Lofstrom T, Kruse A. A Conjunctival Response to Silicone Hydrogel Lens Wear: A new finding reveals how silicone hydrogel lenses may affect the conjunctival epithelium. *Contact Lens Spectrum* 2005; 20(9):42-45.
- Santodomingo-Rubido J, Wolffsohn J, Gilmartin B. Conjunctival epithelial flaps with 18 months of silicone hydrogel contact lens wear. *Eye Contact Lens* 2008; 34:35-38.
- Holden BA, Stephenson A, Stretton S, et al. Superior epithelial arcuate lesions with soft contact lens wear. *Optom Vis Sci* 2001; 78:9-12.
- Pritchard N, Jones L, Dumbleton K, et al. Epithelial inclusions in association with mucin ball development in high-oxygen permeability hydrogel lenses. *Optom Vis Sci* 2000; 77:68-72.
- Tan J, Keay L, Jalbert I, et al. Mucin balls with wear of conventional and silicone hydrogel contact lenses. *Optom Vis Sci* 2003; 80:291-297.
- Szczotka-Flynn LB. Unintended Ortho-k Effects From Silicone Hydrogel Lenses. *Contact Lens Spectrum* 2004; 19(8):19-20.
- Pritchard N, Fonn D, Brazeau D. Discontinuation of contact lens wear: a survey. *Int Contact Lens Clin* 1999; 26:157-162.
- Young G, Veys J, Pritchard N, et al. A multi-centre study of lapsed contact lens wearers. *Ophthalmic Physiol Opt* 2002; 22:516-527.
- Dillehay SM, Miller MB. Performance of Lotrafilcon B silicone hydrogel contact lenses in experienced low-Dk/t daily lens wearers. *Eye Contact Lens* 2007; 33:272-277.
- Dumbleton K, Keir N, Moezzi A, et al. Objective and subjective responses in patients refitted to daily-wear silicone hydrogel contact lenses. *Optom Vis Sci* 2006; 83:758-768.
- Long B, McNally J. The clinical performance of a silicone hydrogel lens for daily wear in an Asian population. *Eye Contact Lens* 2006; 32:65-71.
- Riley C, Young G, Chalmers R. Prevalence of ocular surface symptoms, signs, and uncomfortable hours of wear in contact lens wearers: the effect of refitting with daily-wear silicone hydrogel lenses (senofilcon a). *Eye Contact Lens* 2006; 32:281-286.
- Millodot M, O'Leary DJ. Effect of oxygen deprivation on corneal sensitivity. *Acta Ophthalmol (Copenh)* 1980; 58:434-439.
- Brennan NA, Coles M-LC, Connor HR, et al. Silicone hydrogel contact lens comfort. *Optom Vis Sci* 2007; 84:e-abstract 070037.
- Brennan NA, Coles M-LC. Deposits and symptomatology with soft contact lens wear. *Int Contact Lens Clin* 2000; 27:75-100.
- Ross G, Nasso M, Franklin V, et al. Silicone hydrogels: Trends in products and properties. *Ophthalmic Res* 2005; 2005:27&162.
- Dumbleton KA, Woods CA, Jones LW, et al. Comfort and adaptation to silicone hydrogel lenses for daily wear. *Eye Contact Lens* 2008; 34:215-223.
- Brennan NA, Coles M-LC, Connor HRM, et al. Silicone hydrogel contact lens comfort. *Optom Vis Sci* 2007; 84:E-abstract 070037.